

25th of February 2021

**Confirmation of measurement accuracy**

This is to confirm that the model-type breathalyzers

ACE AL7000

which are determined by

ACE Handels- & Entwicklungs GmbH  
Staufenstr. 1  
83395 Freilassing  
Germany

are always in optimum operating condition and calibrated at the time of delivery. The measuring tolerance of the semiconductor sensor integrated in the ACE AL7000 is +/- 0.10 per mil.

At the time when an ACE AL7000 is delivered, it is ensured that - calculated from the date of shipment of the goods - the next calibration according to the operating instructions of the ACE AL7000 is due either after 300 tests or after 6 (six) to a maximum of 12 (twelve) months - depending on which event occurs earlier.

With kind regards

  
**ACE**  
Instruments  
ACE Handels- & Entwicklungs GmbH  
Staufenstraße 1 · 83395 Freilassing / Germany  
Info@ace-technik.com · www.ace-technik.com  
T+49(0)8654 / 588 38 0 · F+49(0)8654 / 588 38 98

i.A. **Dr. Michael Schuerz**

Head of Alcohol and Drug Measurement Technology (Sales, Marketing & PR)